PLU°LINE
PLURADENT

PLURADES VIRUZID

viruzid

REF 242877

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REF 242877

EN Instrument disinfectant concentrate

Type of application: Product for the manual and terminal disinfection of non-critical and semi-critical medical and surgical instruments in immersion bath procedures. Not suitable for the terminal disinfection of critical medical products!

Notes on application: Dilute the concentrate with 10°C to 30°C warm water (at least drinking water quality) (use dispensing bottles or dispensing equipment). In the event of very dirty conditions, first clean the instruments mechanically. Then immerse the opened instruments in the solution. All surfaces and cavities to be disinfected have to be covered in full by the non-protein attaching solution. Depending on the degree of contamination, the working solution may need to be replaced (RKI recommendation: at least once per working day). After disinfection the instruments should be rinsed, ideally with demineralised water, and dried thoroughly. For critical instruments, sterilisation is recommended in addition as per RKI. Do not mix with cleaning agents or other disinfectants. The preparation guidelines of the instrument manufacturer as well as the preparation guidelines in the RKI/RKINKO recommendation are to be observed! Material changes to polycarbonate and polyamide in the event of long-term use cannot be ruled out, although they can be largely avoided by rinsing sufficiently after each application. When using new materials, the tolerance is to be checked in advance. The working solution can also be disposed of in the sewage water system. To ensure effectiveness, the dosage quantities and application times must be observed as closely as possible and not exceeded. Overdoses should also be avoided for reasons of environmental protection and material tolerance.

Scope of efficacy: Bactericidal (incl. MRSA) as per EN 13727, EN 14561, VAH; levurocidal (candida albicans as per EN 13624, EN 14562, VAH); fungicidal (A. brasiliensis) as per EN 13624, EN 14562; tuberculocidal as per EN 14348, EN 14563; limited virucidal effect (enveloped viruses such as HBV, HIV and HCV) as per DVV/RKI; virucidal (effective against non-enveloped viruses such as norovirus, adenovirus and polio viruses) as per EN 14476.

Application times: Bacteria and Candida albicans (EN 13727, EN 14561, 3624, EN 14562 and VAH), dirty conditions: 2.0 % 15 min, 1.0 % 60 min. TB (M. terre) according to EN 14348, EN 14563 (clean conditions): 2.0 % 15 min, 1.0 % 60 min. Enveloped viruses (e.g. HBV, HCV, HIV) according to RKI/DVV, dirty conditions: 2.0 % 15 min, 1.0 % 60 min. Noro viruses according to

Medizinprodukt/Medical device/
Dispositif médical/Dispositivi medici/
Urządzenie medyczne/Medicinsk indretning/
Medisinsk utstyr/Medicinteknsk produkt

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5Le

CE 0124

UN 1903

UN 1903

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5Le

EN14476 (dirty conditions): 3.0 % 15 min. Adeno viruses according to EN 14476 (clean conditions): 1.5 % 60 min, 3.0 % 15 min. Polyoma SV 40 viruses according to RKI/DVV (dirty conditions): 3.0 % 15 min. Non-enveloped viruses (e.g. Adeno-, Noro-, Polio viruses) according to EN 14476 (clean conditions): 3.0 % 30 min. Fungi (A. brasiliensis) according to EN 13624, EN 14562 (clean conditions): 3.0 % 60 min. Clostridium difficile according to EN 13704 (clean conditions): 3.0 % 60 min.

Safety guidelines: H302 Harmful if swallowed. H314 Causes severe skin burns and eye damage. H410 Very toxic to aquatic life with long lasting effects. Avoid release to the environment. Wear protective gloves / eye protection. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Dispose of contents/container in accordance with local/regional/national/international regulations. Contains diethylmethylammonium chloride, 1-aminopropan-2-ol, N-(3-aminopropyl)-N-decylpropane-1,3-diamine.

Ingredients: Directive 648/2004/EC: Contains 15-15 % non-ionic surfactants, 5-15 % cationic surfactants, <5 % EDTA salts, fragrance.

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PRISMAN GmbH, Otto-Hahn-Ring 6-18, 64653 Lorsch, +49 6251 866980-0, www.prisman.de